



CAMBRIDGE
RESEARCH BIOCHEMICALS

PROTEIN IMMUNISATION SERVICE


(product code CA-14-0050)

Custom polyclonal antibody service

crb provides a simple and cost-effective procedure for producing a rapid polyclonal antibody using a protein

Comprising:

- Y Protein antigen provided by customer (1-2mg)
- Y Guidance on preferred buffer and concentration
- Y Pre-screening of 5 Rabbits pre-immune sera
- Y Customer validation of best 2 Rabbits determined
- Y Two rabbit optimised immunisation program
- Y Free storage for up to 6 months



Pre-screening:

Price £350

- | | |
|----------|--|
| ➤ Day 0 | Receipt of customer own protein and signed protein immunisation guide
5 pre-immune bleeds available for screening |
| ➤ Day 10 | 2 optimum rabbits chosen |

Immunisation Schedule:

Price £960

- | | |
|----------|---|
| ➤ Week 0 | Start of immunisation programme (9 weeks)
Immunisation 1 |
| ➤ Week 2 | Immunisation 2 |
| ➤ Week 3 | Bleed 1 |
| ➤ Week 4 | Immunisation 3 |
| ➤ Week 5 | Bleed 2 |
| ➤ Week 6 | Immunisation 4 |
| ➤ Week 7 | Production bleed 3 |
| ➤ Week 8 | Immunisation 5 |
| ➤ Week 9 | Final harvest bleed |

Additional Services:

- | | |
|---------------|---|
| Code: | CA-14-0020 |
| Service: | Elisa analysis of all bleeds |
| Price: | £80 per bleed |
| Code: | CA-14-0030 |
| Service: | Affinity or protein A purification of final harvest bleed |
| Price: | £590 |
| Code: | CA-14-0040 |
| Service: | Charge for each additional rabbit |
| Price: | £480 |

Pre-screening: £350 + Immunisation: £960
Total price: £1310

PROTEIN IMMUNISATION GUIDE

Protein

Antigen name

How is the protein derived?

HINT: Host, expression system, soluble/inclusion body, refolded protein, purification method (e.g. 6-His, GST or MBP Tag). If you have sub-contracted the expression of the protein, please provide supplier's details and datasheet.

How is the protein to be supplied?

HINT: **In solution** If so, what is the buffer?
Lyophilised If so, which buffer should we re-suspend it in?
PAGE slice

Total antigen to be provided (mg)

As a guideline, we need 1-2mg of protein for a 2-rabbit immunisation programme.

Buffer

In general, most buffers and additives are suitable for antibody production apart from detergents, protein solubilising and refolding agents such as urea and guanidinium hydrochloride and organic solvents such as DMSO.

Immunisation

Which species do you require and how many?

species

number

HINT: We would typically recommend a minimum of 2 rabbits. Other species include rat, mouse, goat, sheep, chicken.

Screening for Antibody

For ELISA analysis, we need an additional 1mg of antigen to coat immunoassay plates. It would be better to provide an untagged protein or differently tagged protein so that we can measure the immune response specifically to the target protein (and not the tag). We can assay tagged protein, however, and additionally perform ELISA screening for antibodies directed against the tag itself as a comparative control.

Animal Welfare Assurance

I CERTIFY THAT THE ANTIGEN I AM SENDING TO CAMBRIDGE RESEARCH BIOCHEMICALS LIMITED FOR THE PURPOSES OF RAISING A POLYCLONAL ANTIBODY IS FOR RESEARCH PURPOSES ONLY, AND IS **NOT DANGEROUS TO PEOPLE OR ANIMALS. FURTHERMORE, THE SUPPLIED ANTIGEN IS NOT A LIVE BACTERIA/VIRUS, ANIMAL/HUMAN PATHOGEN OR TOXIC/HAZARDOUS TO ANIMALS OR PERSONNEL**

Signature		Print Name	
Email		Date	
Company / Insititute Adress			

Important Note

Please provide a datasheet with details of protein name, molecular weight (kDa), amount (mg), buffer, concentration (mg/ml) and volume (ml) per tube. Please label tube(s) clearly with volume & concentration or weight.

Fax this completed form to +44 (0)1642 567181
OR scan & email to immunology@crbdiscovery.com

Immunology Department

Cambridge Research Biochemicals Limited,
17-18 Belasis Court, Belasis Hall Technology Park,
Billingham, Cleveland TS23 4AZ United Kingdom

T +44 (0) 1642 567180 F + 44 (0) 1642 567181
E enquiries@crbdiscovery.com
www.crbdiscovery.com

